









Biohazardous Waste Do's and Don'ts

| Do | Don't |
|---|--|
|  <p>Dispose of blunt tip pipettes in a biomedical waste container.</p> |  <p>Don't pre-treat disposable pipettes in a graduated cylinder.</p> |
|  <p>Use properly labeled biomedical waste containers with the biohazard symbol or the words "Biohazardous Waste".</p> |  <p>Don't use non-labeled containers.</p> |
|  <p>Tie off biomedical bag when $\frac{3}{4}$ full.</p> |  <p>Don't overfill container.</p> |
|  <p>Keep lid closed except when adding waste.</p> |  <p>Don't keep lid open.</p> |

Do not dispose of biohazardous waste in the regular trash.